



DATA SHEET

Brilliant Violet 785 anti-mouse CD4 Antibody

Code Number#21-100453

Pack Size: 50 µg

| | |
|--------------------------------|--|
| Clone: | GK1.5 |
| Isotype: | Rat IgG2b, κ |
| Reactivity: | Mouse |
| Antibody Type: | Monoclonal |
| Host Species: | Rat |
| Immunogen: | Mouse CTL clone V4 |
| Formulation: | Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA). |
| Preparation: | The antibody was purified by affinity chromatography and conjugated with Brilliant Violet 785 under optimal conditions. |
| Concentration: | 0.2 mg/ml |
| Application: | FC - Quality tested |
| Storage & Handling: | The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze. |

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.