

DATA SHEET**ALK-3 / BMPR1A Antibody, Rabbit PAb, Antigen Affinity Purified****Code Number#21-100743-T08**

Immunogen:	A synthetic peptide corresponding to the N-terminus of the Human ALK-3 / BMPR1A
Preparation:	Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human ALK-3 / BMPR1A, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human ALK-3 / BMPR1A
Formulation:	0.2 µm filtered solution in PBS
Applications:	IHC-P
Storage:	This antibody can be stored at 2°C -8°C

For in-vitro diagnostic Use. Not For Human Or Animal Consumption.