

**DATA SHEET****GLP1 Antibody****Code Number#21-NBP2-23561BTL****Pack Size: 0.1 ml**

<b>Clonality:</b>	Monoclonal
<b>Clone:</b>	11E2
<b>Specificity/Sensitivity:</b>	Reacts with all forms of GLP-1, including precursor and GLP-1(9-37) /GLP-1(9- 36 amide) metabolite
<b>Immunogen:</b>	Synthetic GLP-1(7-36) amide coupled to carrier and adsorbed onto aluminum hydroxide gel
<b>Preservative:</b>	0.05% Sodium Azide
<b>Isotype:</b>	IgG2a Kappa
<b>Conjugate:</b>	Biotin
<b>Purity:</b>	Protein A or G purified
<b>Buffer:</b>	PBS
<b>Gene Symbol:</b>	GCG
<b>Storage:</b>	Store at 4C in the dark.

**Application Note:** Optimal dilution of this antibody should be experimentally determined.**Note: For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.**