

DATA SHEET**Anti-hPD-L1-hIgG1, purified antibody****Code Number#21-hpdl1-mab1-01****Pack Size: 100 mg**

Target:	Human programmed cell death ligand 1 (PD-L1)
Species reactivity:	Reacts with human and mouse PD-L1
Isotype:	Human IgG1
Light chain type:	kappa
Clonality:	Monoclonal
Source:	CHO cells
Purity:	Protein G purified
Formulation:	0.2 mm filtered solution in a sodium phosphate buffer with glycine, saccharose and stabilizing agents
Applications:	Flow cytometry
Antibody resuspension:	Add 1 ml of sterile water to obtain a concentration of 0.1 mg/ml.
Storage:	Lyophilized product is stable for 1 year when stored at -20 °C.

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.

