

**DATASHEET****Goat anti-Mouse IgG Antibody (HRP)**

Code Number#336668-08KT

<b>Target:</b>	IgG
<b>Reactivity:</b>	Mouse
<b>Host:</b>	Goat
<b>Clonality:</b>	Polyclonal
<b>Conjugate:</b>	HRP
<b>Application:</b>	Immunoassay (IA), Immunohistochemistry (IHC), ELISA, Western Blotting (WB)
<b>Immunogen:</b>	Goat antiserum was obtained from animals of US origin
<b>Characteristics:</b>	May contain small amounts of intact IgG
<b>Purity:</b>	> 90 % based on SDS-PAGE
<b>Applications:</b>	This conjugate is suitable for all immunoassay applications. The optimal working dilution should be determined by the investigator. Suggested starting dilution: 1:500-1:5,000 for Immunohistochemistry 1:200-1:5,000 for ELISA/Western blot
<b>Storage:</b>	4 °C/-20 °C

**For In-Vitro Diagnostic Use only. Not for human or animal consumption.**