

DATA SHEET**Myelin Oligodendrocyte Glycoprotein 35-55 (MOG 35-55)****Code Number#BY-3038001-90****Pack Size:** 25 mg/vial**Description:** Purified Myelin Oligodendrocyte Glycoprotein (MOG) 35-55.**Form:** Lyophilized**Purity:** $\geq 96\%$ (HPLC)**Molecular Weight:** 2581.96**Counterions Present:** TFA salt**Application:** Multiple Sclerosis: MOG 35-55 induced experimental autoimmune encephalomyelitis (EAE)**Storage:** Shipping at ambient temperature. Storage at -20°C .

Note: For In-Vitro Diagnostic Use only. Not for human or animal consumption.