

**DATA SHEET****Anti-Urokinase, Polyclonal Rabbit**

Code Number#CAT91-605

|                         |   |
|-------------------------|---|
| <b>Isotype:</b>         | Mostly IgG.   |
| <b>Titre:</b>           | 1:100 to 1:10,000 depending on the application.   |
| <b>Antigen:</b>         | Two chain human urinary type high molecular weight plasminogen activator isolated from pooled urine.  |
| <b>Cross-reactions:</b> | Reacts only with human u-PA isolated from kidney cell culture tissue.<br>No reactivity with other human proteins. Identity established by double diffusion with u-PA isolated from human kidney tissue and E. coli- derived recombinant u-PA. No cross reactivity with human tissue type plasminogen activator. |
| <b>Characteristics:</b> | Inhibits urokinase activity and recognises u-PA-PAI-1 and -2 complexes.   |
| <b>Presentation:</b>    | 1ml unpurified serum and no preservatives.  |
| <b>Applications:</b>    | Activity inhibition characterisation (especially coupled with fibrin autography) and RIA. Immunoprecipitation, immunocytochemistry using frozen or paraffin-fixed sections, Western Immunoblotting, ELISA, fluorescence activated cell sorting.   |
| <b>Storage:</b>         | +4°C for 1 month provided preservative is added For extended storage, freeze in working aliquots at -20°C. Avoid repeated freezing and thawing.   |

**For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.**