



## DATA SHEET

**Cell-Surface Vimentin (CSV) monoclonal antibody, clone 84-1**

**Code Number#H00007431-M08**

**Description:** Mouse monoclonal antibody recognizes human cell-surface vimentin (CSV).

**Clone Name:** 84-1

**Immunogen:** Human recombinant vimentin

**Host:** Mouse

**Reactivity:** Human

**Applications:** Flow Cyt, IF

**Form:** Liquid

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 7431

**Gene Symbol:** VIM

**Gene Alias:** FLJ36605

**For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.**