

DATA SHEET**Anti-Streptokinase (Polyclonal Rabbit)****Code Number#KT-CA-1012b-029**

- Immunogen:** Highly purified Streptococcus haemolyticus.
- Titre:** The titre value is at least 1/32 by double immunodiffusion or 1/1000 by ELISA.
- Presentation:** Unpurified serum.
- Cross Reactivity:** These antibodies have single antigen reactivity with minimal cross reactivity with other human serum proteins.
- Applications:** Suitable for a variety of techniques including:
- a) Detection of target antigens by: Double Immunodiffusion / Immunoelectrophoresis / 2-D IEP
 - b) Quantitative Assay by: ELISA / Laurel Rocket IEP
 - c) Immunoinhibition
- Storage:** Store at +4^oC in the short term. Otherwise aliquot and store at -20^oC.

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.