

DATA SHEET**POOLED HUMAN NORMAL PLASMA****Code Number#IPLA-K2**

- Description:** Normal Human Plasma is processed from single donors via apheresis plasma units. The material undergoes viral testing, and the final product is aliquoted to customer specifications.
- Viral Testing:** Each unit is tested and found negative for HBsAg, HCV, HIV-1, HIV-2, HIV-1Ag or HIV 1-NAT, ALT, and syphilis by FDA-approved methods.
- Expiration:** Based upon industry standards, plasma is stable for 3 years from the date of manufacture.
Please see product label for specific date of expiration. Pooled Normal Human Plasma should be stored at -20° to -80° Celsius.
- Handling:** This material is sold for For In-Vitro Diagnostic Use only. Not for human or animal consumption. It is the responsibility of the user to undertake sufficient verification and testing to determine the suitability of each product's application.
- The statements herein are offered for informational purposes only, and are intended to be used solely for your consideration, investigation, and verification.