

DATA SHEET**Rubeola Nucleoprotein Antigen, Recombinant**

Code Number#K81-EV9298

Viral Strain:	Edmonston
Source:	Grown in <i>Sf-9</i> Cells Product manufactured in the United States using only domestic animal materials processed in U.S. facilities.
Format:	Purified Liquid
Protein Concentration:	BCA Assay:
DOM:	Acceptable Range: 0.4 – 0.6 mg/mL Testing Results: 0.49 mg/mL Testing Results: Pass
Identity:	SDS-PAGE
Western Blot:	Pass
Buffer:	50 mM HEPES, pH 7.4 – 7.6
Preservative:	None
Applications:	ELISA
Storage:	Store at $\leq -65^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles

FOR IN-VITRO DIAGNOSTIC USE. NOT FOR HUMAN OR ANIMAL CONSUMPTION