



DATA SHEET

Antibodies Hepatitis E Virus Code

Number#KLS-C67675-025

Description:	Hepatitis E Virus antibody LS-C67675 is an unconjugated mouse monoclonal antibody to hepatitis e virus Hepatitis E Virus (HEV) (aa 434-457). Validated for ELISA, IF and WB. Cited in 1 publication.
Target:	HEV / Hepatitis E Virus
Synonyms:	HEV
Host:	Mouse
Reactivity:	Hepatitis E Virus(tested or 100% immunogen sequence identity)
Clonality:	IgG2b Monoclonal
Conjugations:	Unconjugated. Also available conjugated with AP, APC, Biotin, FITC, HRP, PE.
Purification:	Purified
Modifications:	Unmodified
Immunogen:	Hepatitis E recombinant protein encoding amino acids 267-660 of ORF2.1.
Specificity:	ORF2.1, encoded region to capsid protein, the epitope is mapped to aa 434-457.
Presentation:	0.02M PBS, pH 7.6, 0.09% sodium azide.
Storage:	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.