

**DATA SHEET****Pichia Pastoris Antibody, polyclonal****Code Number#KT-20R-PG008**

<b>Species:</b>	Yeast
<b>Expression System:</b>	Rabbit / IgG
<b>Class:</b>	Polyclonal
<b>Type:</b>	Antibody Immunogen Whole cells of Pichia pastoris.
<b>Conjugate:</b>	Unconjugated
<b>Form:</b>	Liquid
<b>Concentration:</b>	4-5 mg/mL
<b>Purification:</b>	Purified
<b>Storage Buffer:</b>	PBS, pH 7.2, Contains 0.1% sodium azide
<b>Applications:</b>	Immunofluorescence (IF), Western Blot (WB)
<b>Storage Conditions:</b>	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C

**For In-Vitro Diagnostic Use only. Not for human or animal consumption.**