



DATA SHEET

Anti -AR-V7 Monoclonal Mouse Antibody

Code Number#KT-AG10008

Androgen Receptor: Variant 7

Isotype: IgG 2a

Cross-Reactivity: Human

Pack Size: 100 µg (1mg/ml) purified antibody in PBS, pH 7.4, or Lyophilized.

Applications: Western Blot 1–2 µg/ml I HC1–3 µg/ml. However, the optimal conditions should be determined individually.

Storage: Store at -20°C. Avoid repeated freezing and thawing.

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption