



DATA SHEET

GLP1 Antibody (5B10) - BSA Free

Code Number: KT-NBP1-05189

Unit Size: 100 ug

Concentration: 1.0 mg/ml

Storage: Store at 4C. Do not freeze.

Clonality: Monoclonal

Clone: 5B10

Preservative: 0.05% Sodium Azide

Isotype: IgG1 Kappa

Purity: Protein A or G purified

Buffer: 10mM Phosphate (pH 7.4) and 0.5M NaCl

Host: Mouse

Reactivity Notes: The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

Specificity/Sensitivity: Reacts with all forms of GLP-1, including precursor and GLP-1(9-37) /GLP-1(9- 36 amide) metabolite.

Immunogen: Synthetic GLP-1(7-36)amide coupled to carrier and adsorbed onto aluminum hydroxide gel.

Applications: Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-

Note: This product is For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.

Primary Antibodies are guaranteed for 1 year from date of receipt.