



## Certificate of Analysis

### Pooled Human Male AB Serum - Heat Inactivated

Cat No:	KT-SERAB-HI
QTY:	SAMPLE
SIZE:	100ml / 1000ml
Store:	-20° Celsius
Expiration:	SAMPLE
Glucose	96 mg/dl
Sodium	140 meq/l
Potassium	3.7 meq/l
Chloride	120 meq/l
Calcium	8.3 mg/dl
Phosphorus	2.6 mg/dl
Total Protein	5.3 gm/dl
Albumin	3.3 gm/dl
Total Bilirubin	0.3 mg/dl
Cholesterol	130 mg/dl
Triglycerides	72 mg/dl
Hemoglobin	6.3 mg/dl
pH	7.6
Osmolarity	286 mOsm
Sterility	Sterile
Mycoplasma	Negative
Endotoxin	0.8 EU/ml Growth Promotion Passes
Physical Appearance	Clear, straw yellow solution Filtration 0.1µm sterile filter

Each donor unit used in the preparation of this material was tested and found negative for HBsAG, HIV-1/2, anti-HCV, HCV-Nat, HIV-1 NAT, and Syphilis per FDA Regulations. This product is collected and processed from male donors, blood type AB, in FDA-licensed commercial donor establishments and FDA-registered facilities located in the USA.

This material is sold for in-vitro diagnostics use. Not for human or animal consumption. This material is not suitable for human use. It is the responsibility of the user to undertake sufficient verification and testing to determine the suitability of each product's application. The statements herein are offered for informational purposes only and are intended to be used solely for your consideration, investigation, and verification.

Caution to Users: The information expressed in this product data sheet represents best efforts to provide accurate and reproducible data about the material. However, test methods and techniques for biologicals may produce variable results; therefore, the user should evaluate the material for suitability in a given system prior to purchase. Lot specific samples are available at the customer's request and users are encouraged to evaluate the sample prior to purchase.

**For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption**