



DATA SHEET

ChromPure Mouse IgG, whole molecule

Code Number: #KT91015-000-003

Lot Number: 151748

Physical State: Sterile-filtered liquid

Size: 5.0mg

IgG Concentration: 5.7 mg/ml

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH7.6

Preservative: None

Storage: Store at 2-8°C under sterile conditions. Prepare working dilution on day of use.

Expiration Date: one year from date of receipt.

Purity: Based on immunoelectrophoresis at an antigen concentration of 20 mg/ml, the pattern of precipitation against anti-mouse whole serum is the same as that against anti-mouse IgG, Fc fragment specific.

Note: For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.