



## DATA SHEET

### Goat Anti-Human IgG, Fc Fragment Specific

**Code Number#**KT91109-001-008

**Lot Number:** 150134

**Physical State:** Freeze-dried solid

**Size:** 2.0 ml

**Protein Concentration:** 33.7 mg/ml

**Buffer:** 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

**Preservative:** 0.05% Sodium Azide

**Storage and Rehydration:** Store freeze-dried solid at 2-8°C. Rehydrate with 2.0 ml dH<sub>2</sub>O and centrifuge if not clear.

**Expiration Date:** one year from date of rehydration. The expiration date maybe extended if test results are acceptable for the intended use.

**Purity:** The antibody was delipidated and ammonium sulfate fractionated from whole antiserum.

**Note:** For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.