



DATA SHEET

IL-2 Monoclonal Antibody (JES6-5H4), APC

Code Number#KT9117-7021-82

Size:	100 µg
Species:	Mouse
Class:	Monoclonal
Type:	Antibody
Clone:	JES6-5H4
Conjugate:	APC
Form:	Liquid
Concentration:	0.2 mg/mL
Purification:	Affinity chromatography
Storage Buffer:	PBS, pH 7.2, with 0.1% gelatin Contains 0.09% sodium azide
Storage Conditions:	4° C, store in dark, DO NOT FREEZE!
RRID	AB_469490
Applications:	Flow Cytometry (Flow), ELISA (ELISA)

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption