



DATA SHEET

Human Parathyroid Hormone (PTH N- terminal) Monoclonal Antibody

Cat Number: KTA66030

Source: Tissue Culture Supernatant
Immunogen: Synthetic PTH (aa1-34)
Purification: Protein A affinity chromatography
Buffer: Phosphate buffered saline, pH 7.2
Preservative: 0.1% sodium azide
Subclass: IgG1
Storage: 2-8 degrees short term, -20 degrees long term
Specificity: This antibody reacts with 1-34 PTH

This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption