



## DATASHEET

### Anti-Exendin 4 antibody

Code Number#KTAB106796-02

<b>Host species:</b>	Rabbit
<b>Specificity:</b>	Exendin 4 is a 39 amino acid peptide found in venom from the Gila monster Heloderma suspectum.
<b>Species reactivity:</b>	Reacts with: Heloderma suspectum
<b>Immunogen:</b>	Synthetic peptide: HEGGTFTSDL SKQMEEEEAVR LFIEWLKNGG PSSGAPPPS conjugated to KLH via carboxyl group.
<b>Form:</b>	Liquid
<b>Storage instructions:</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Storage buffer pH: 7.50, Preservative: 0.02% Sodium azide Constituents: PBS, 50% Glycerol
<b>Purity:</b>	Protein G purified
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG

**Note: All products are "For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption."**