

DATA SHEET**RABBIT ANTI-VASOPRESSIN POLYCLONAL ANTIBODY****Code Number**#KTAB1565**Pack Size:** 100 µL**Specificity:** Vasopressin. Less than 1% cross-reactivity with oxytocin.**Immunogen:** Arginine vasopressin (Sigma) conjugated to thyroglobulin.**Applications:** Immunohistochemistry, Immunoblotting, RIA**Sensitivity:** < 1 pg/mL.**Species Reactivity:** Human, rat, mouse, sheep and rabbit.**Format:** Neat rabbit antisera.**Presentation:** Liquid, with 0.1% sodium azide.**Storage/Handling:** Store at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.**For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.**