

DATASHEET**Anti-Exendin 4 antibody****Code Number: KTAB2340713**

Host species:	Mouse
Species reactivity:	Reacts with: Recombinant fragment
Immunogen:	Carrier coupled exendin 4.
Epitope:	The epitope is located in the region of amino acids 9 to 39 of exendin 4.
Form:	Liquid
Storage instructions:	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer:	pH: 7.40, Preservative: 0.097% Sodium azide, Constituents: 0.0268% PBS, 2.9% Sodium chloride
Purity:	Protein A purified
Clonality:	Monoclonal
Clone number:	20
Isotype:	IgG1

Note: All products are "For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption."