


**DATA SHEET****Calprotectin, protein**

Code Number#KTATC04-02

<b>Description:</b>	Heterotetramer containing two 8 kDa and two 14 kDa subunits
<b>Source material:</b>	Human neutrophils
<b>Purity:</b>	>90% by SDS-polyacrylamide electrophoresis
<b>Purification:</b>	purified using conventional chromatography procedures
<b>Clinical Applications:</b>	Fecal calprotectin determination has been shown to be useful in the diagnosis and management of inflammatory bowel disease



SDS-polyacrylamide electrophoresis of ATC04 Calprotectin antigen showing two subunits of m.wt. 8 and 14 kDa

**For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption**