



DATA SHEET

Anti-Bevacizumab monoclonal antibody, clone 2C8. Mouse monoclonal antibody raised against F(ab)2 fragment of Bevacizumab.

Code: KTDS-MAB11128

Clone Name: 2C8-13

Immunogen: F(ab)2 fragment of Bevacizumab.

Host: Mouse

Applications: Blocking, ELISA

Specificity: This antibody blocks the Bevacizumab-VEGF interaction in ELISA.

Form: Liquid

Purification: Protein G purification

Isotype: IgG1, kappa

Recommend Usage: ELISA (2 ug/mL for coating) The optimal working dilution should be determined by the end user.

Storage Buffer: In TBS, pH 8.0 (0.09% sodium azide)

Storage Instruction: Store at 4°C for at least 12 months. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption