



DATA SHEET

Mumps Virus Antibody

Code Number#KTHM050

Product Type: Monoclonal

Pack Size: 1 mg

Format: Protein A Purified Monoclonal Antibody

Host/Source: Mouse

Purity: >90%

Specificity: Nucleoprotein. No cross-reactivity with Measles Virus, RSV, Parainfluenza 1, 2, 3

Class/Isotype: IgG2b

Supplied In: 10 mM Phosphate Buffered Saline, pH 7.2

Preservative: 0.1% Sodium Azide

Applications: ELISA, Lateral Flow, IFA

Storage: Short Term at 2-8°C. For Long Term storage at -20°C. Avoid repeated freezing and thawing.

Note: This product is intended for in-vitro diagnostic use. not for human or animal consumption.