

DATA SHEET

Mumps Virus Antibody

Code Number#KTHM400

Product Type:	Monoclonal
Pack Size:	1 mg
Format:	Purified Monoclonal Antibody
Host/Source:	Mouse
Purity:	>95% by SDS-PAGE
Immunogen/Antigen:	Purified Mumps Virus
Class/Isotype:	IgG
Supplied In:	10 mM Phosphate Buffered Saline, pH 7.0
Preservative:	0.02% Sodium Azide
Applications:	ELISA, Lateral Flow
Storage:	Short Term at 2-8°C. For Long Term storage at -20°C. Avoid repeated freezing and thawing.

Note: This product is intended for in-vitro diagnostic use. not for human or animal consumption.