

DATA SHEET**Histone H3K4me3 Polyclonal Antibody****Code Number#KTHW0031**

Host Species:	Rabbit
Clonality:	Polyclonal
Purification:	Antibodies were produced by immunizing rabbits and were purified by antigen affinity-chromatography.
Applications:	WB, IHC, IF, IP, ChIP
Species Reactivity:	Hu, Ms, Rt
Immunogen:	Type Peptide
Immunogen Description:	A synthetic methylated peptide corresponding to residues surrounding K4 of human histone H3
Target Name:	Histone H3
Modification:	Methyl
Concentration:	1.0mg/ml
Formulation Buffer:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage:	Store at -20°C or -80°C. Avoid freeze / thaw cycles.
Modification:	Methyl
Concentration:	1.0mg/ml
Formulation Buffer:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage:	Store at -20°C or -80°C. Avoid freeze / thaw cycles.

Note: For In-Vitro Diagnostic Use. Not For Human Or Animal Consumption.