



## DATA SHEET

### Hepatitis A Virus (HAV) Antibody

Code Number#KTJ128

<b>Product Type:</b>	Polyclonal
<b>Pack Size:</b>	500 µl
<b>Format:</b>	Purified Polyclonal Antibody, Ig Fraction
<b>Host/Source:</b>	Goat
<b>Purity:</b>	>95%
<b>Immunogen/Antigen:</b>	Purified HAV HM175
<b>Specificity:</b>	Virions
<b>Supplied In:</b>	10 mM Phosphate Buffered Saline pH 7.2
<b>Preservative:</b>	0.1% Sodium Azide
<b>Applications:</b>	ELISA, IFA, WB
<b>Storage:</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Note: This product is intended for in-vitro diagnostic use. Not for human or animal consumption.**