

DATA SHEET**GLP1 Antibody (11E2) [Biotin]****Code Number#KTNBP2-23561B47**

Pack Size:	0.1 ml
Clonality:	Monoclonal
Clone:	11E2
Preservative:	0.05% Sodium Azide
Isotype:	IgG2a Kappa
Conjugate:	Biotin
Purity:	Protein A or G purified
Buffer:	PBS
Host:	Mouse
Gene ID:	2641
Gene Symbol:	GCG
Reactivity Notes:	Mammal
Specificity/Sensitivity:	Reacts with all forms of GLP-1, including precursor and GLP-1(9-37) /GLP-1(9- 36amide) metabolite
Immunogen:	Synthetic GLP-1(7-36)amide coupled to carrier and adsorbed onto aluminum hydroxide gel
Storage:	Store at 4C in the dark.
Application:	Western Blot, ELISA

Note: For in-vitro diagnostic use only. Not for human or animal consumption.